




ResMed masks with magnets

Note: product availability may differ in each country.



ResMed masks with magnets have magnetic clips to facilitate a simple and easy way for patients to attach and detach the headgear to the mask frame.

Awaken your best.

www.resmed.com/magnetupdate

	Location of model name	Location of magnets
AirFit™ F30 Full face mask		
AirFit™ F30i		
AirFit™ F20 Full face mask AirTouch™ F20 Full face mask AirFit™ F20 Full face mask for Her AirTouch™ F20 Full face mask for Her	 	
AirFit™ N20 Nasal mask AirTouch™ N20 Nasal mask AirFit™ N20 Nasal mask for Her AirTouch™ N20 Nasal mask for Her	 	
AirFit™ N10 Nasal mask AirFit™ N10 Nasal mask for Her		
AirFit™ F20 NV		